6 in 1 Ultrasonic 40K Cavitation RF & Laser Slimming Unit



40K Cavitation Head

Power supply output : 150v Frequency : 400 Watt Power : 25 Watt

Vacuum + Bipolar RF Head

Power supply output : 66v Frequency : 1 MHz Power : 95 Watt





Face RF Head Power supply output: 66v

Frequency : 450 khz Power : 35 Watt



Body RF Hea

Power supply output Frequency : 550 khz Power : 25 Watt



Laser Pads

Lenght: 635 - 650nm Power: 400 Watt Single output: 50 mw Number of diode: 60 diod:



Energy output	56 x 180 mW
Main power output	100 VA
Safety	On board diagnostics
Classification	Electrical Class 1, Type B applied part
Cooling requirements	Air cooled
Mode of operation	Continuous



6 in 1 Slimming Machine

RF and ultrasonic beauty machines are currently the most popular beauty machines that work on tightening, body-shaping and anti-aging. It's an effective replacement of all the facial and body care programs in beauty salons. Moreover, they are convenient and easy to operate. It solves both skin and body management issues for beauty-lovers.

This equipment requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. RF machines heats deeper skin layers, stimulates collagen cells to reform in dermis, increases blood circulation in hypodermis, thus achieving skin tightening and anti-aging.

Advantages

- 6-in-1 multifunctional beauty machine that works on both of your facial and body care.
- 2. Ultrasonic fat burning replaces and transcends liposuction, bringing good news for treating obesity.
- 3. Various work heads can be replaced in different parts according to different requirements.
- 4. The perfect combination of ultrasound and RF solves the dilemma of losing weight. It helps with weight loss while at the same time maintains firm skin, increases skin elasticity, metabolism, and ultimately better health.
- 5. Adopting the cutting-edge technology of ultrasonic fat burning.
- 6. Suitable for all skin types.
- 7. Painless, non-invasive during the process. Skipping any recovery periods without affecting any work and life plans.
- 8. No consumption, low cost and quick returns.
- 9. Wider treatment range and faster and more visible effects.
- 10. Unevenness, bleed, swelling and stasis will not appear after treatment.

Handles included:

40K Cavitation Head:

Cavitation causes vibration of the liquid parts inside a fat cell, the vibration erupts the fat cell from the inside out, exposing the fat, and eliminating the fat pocket and membrane.

RF (Radio Frequency) Handle:

Radio frequency tones, tightens and stimulates collegen in the skin. Lifting and firming of areas such as the buttocks, thighs, breasts, face and eyebags. This unit has a small, medium and large handle.

Vacuum Handle:

Lifts areas and is used to soften the appearance of cellulite.

Laser Lipo

Cold laser lipo ensures to liquify the free fat waste, into a thin oily liquid, making lymph drainage easier after a cavitation session, used on its own, Lipo laser is effective to reduce a fat cell size, as it creates microscopic injuries to a fat pocket, and in turn drains the fat cell.

N\$19 798.90

Slimming Technology that works