

The main role of the low-energy laser bio-stimulation (Biostimulator), is, through the give appropriate energy to stimulate the biological cells and induce or strengthen several physiological responses, including promoting local blood circulation, regulating cell function, enhancing immune function, promoting cell metabolism and proliferation.

Wavelength of 650nm-660nm red laser just in human eye color of the visible spectrum, so we can see the red light 650nm -660nm can penetrate the organization up to 8-10mm, the effective activation and repair cells, stimulate the synthesis of cell metabolism, for the superficial cells biochemical stimulation and hyperemia. Irradiation meridian points to stimulate the meridian points related organizations, the skin tissue will not cause harm to exempt patients fear of needles, and also have the function to stimulate the meridians in a healthy way.

#### 635nm

Fat reduction, increase cell-energy, Mitochondrien and ATP production, to supply body energy and promote metabolism.

#### 660nm

Stimulate the lymph drainage and reduce the cellulite.

#### 810nm

For all bone-illnesses (Arthrosis, arthritis, carpal tunnel), it stimulates reparative cells and microcirculation.

#### 830nm

The laser in the 830nm band penetrates deeper into adipose tissue, it accelerates the metabolism of cellulite, promotes rapid aging of cellulite, and excretes it from the body through human metabolism.

#### 940nm

Invisible light stimulate the lymph drainage and reduce the cellulite, more effective, penetrates deeper into adipose tissue.

#### 980nm

Same with 660nm, more effective to reduce fat and tighten the skin.

## Specifications

Laser source	Mitsubishi Diode Laser ML101J27
Laser type	Low Intensity Diode Laser
Quantity of lamps	528 pcs
Lamp energy	100 mW - 350 mW (Pulse)
Laser working modes	Continuous and Pulse
Laser pads	12 pcs (8 Large Pads and 4 Small)
Touch screen	8 Inch Touch Screen
Cooling system	Air cooled
Voltage	110V/220V; 50Hz - 60Hz
Safety	Onboard diagnostic
Classification	Electrical Class II, Type B Applied part
Environmental requirement	10°C to 40°C. Non-condensing humidities below 75% RH
Package size	49 x 66 x 59cm
Gross Weight	30 kg

# Portable 6D LASER LIPO Slimming Machine




Product available in black & white



12 Laser Pads Included

Fat-Measure Caliper



Protective Goggles



**635nm, 660nm, 810nm, 830nm, 940nm and 980nm combined in one amazing system.**

**N\$153 998.90**